

Method development for high-throughput sample analysis



The challenge

CatSci was tasked with catalyst screening for an asymmetrical hydrogenation of a target. To quickly explore reaction conditions, high throughput screens (up to 96 experiments per single run) were carried out. The original customer method had a run time of 35 minutes – limiting for high throughput development. A new faster method was required that maintained the required critical resolution between peak pairs.




How?

Column and chromatographic conditions were screened using a design of experiments approach, to optimise method parameters. CatSci's extensive analytical development experience maintained critical pair resolution and delivered methods with dramatically shorter analysis time.



The achievement

A 6-minute UPLC method was successfully developed which allowed the chemist to analyse the catalyst screening results with 24 hours of completing the experiment.



Resolution of critical pairs maintained

Old Process

Achiral UPLC screening method – 35 minute analysis time

Potential analysis time for 96 samples = 56 hours

Amount of solvent A (aqueous buffer) used for a 96 sample run = 495 mL

Amount of solvent B (organic solvent) used for a 96 sample run = 512 mL

New Process

Achiral UPLC screening method – 6 minute analysis time

Potential analysis time for 96 samples = 9.6 hours

Crystallisation from reaction solvent

Higher yields and consistent purity

Facts and Figures

Water Waste
44% less



Solvent Waste
36% less



Analysis time saved
83% less



AMVI*
42% less



*Analytical Method Volume Intensity (Green Metric)



The team: Dr Siân Forsyth, Gareth Tutton & Beth Rees



Find out more: catsci.com